


LIRIA

1000 GPD



DOMESTIC REVERSE OSMOSIS DIRECT FLOW

The RO LIRIA thanks to its 1000 GPD membrane is able to produce a large quantity of pure water in a short time. This makes it ideal for homes with high demand, always ensuring an efficient and continuous supply. Power and purity in every drop.



FAST CONNECTION
Lifts and extracts filter and membrane. Safe, efficient and fast mechanics.



DIRECT FLOW
No tank required.



STATUS INDICATOR LEDES
Information on cartridge life.



HIGH EFFICIENCY
High-performance technology that achieves maximum flow with less rejection.



CARTRIDGE AND MEMBRANE
Original and exclusive



1000 GPD
High performance for high flow rates.



SMART FLUSHING
Different types of automatic flushing (first use, accumulated operation, change of cartridges...).



AUTOMATIC MAINTENANCE NOTIFICATION
Through LEDs and sound



QUICK AND EASY MAINTENANCE



European Certificate of Conformity


LIRIA

1000 GPD


REVERSE OSMOSIS DIRECT FLOW

Filtering system containing:
1.- PAC Pretreatmet Cartridge
2.- ME 1000 GPD Membrane

- Intuitive interface with led prompts for cartridge and membrane.
- Includes tap
- Peak flow rate: 2,62 Lpm
- Membrane 1000 GPD
- Maximum inlet TDS: 1.000 ppm
- High performance diaphragm booster pump
- Inlet Pressure: 1 to 4 Bar
- Inlet temperature: 4 to 38 °C
- Power supply: 24Vdc.
- External electrical adapter: 100-240V 50/60Hz: 24Vdc



EQUIPMENT MEASUREMENTS				
CODE	DIMENSIONS (mm)			
	A	B	C	PESO
RO-8100-02	145	380	385	10 Kg



PACKAGING AND TRANSPORT				
CODE	DIMENSIONS (mm)			LOGISTICS
	A	B	C	
RO-8100-02	578	260	486	Pallet 24 uds.

HIGH RANGE DIRECT FLOW Reverse Osmosis:

- Quick and easy connection
- 1000 GPD: high flow rate
- No tank
- Does not take up space under the sink
- European CE Conformity Certificate



CODE	DESCRIPTION
CA-6100-02	PAC Pretreatmet Cartridge
ME-6100-02	ME 1000 GPD Membrane
RO-8100-02	Liria 1000 GPD Direct Flow Reverse Osmosis